


Cardio Flow

FreedomFlow® Peripheral Guidewire

Designed to tackle lesions both above and below the knee

FreedomFlow®
by Cardio Flow
PERIPHERAL GUIDEWIRE

The FreedomFlow Peripheral Guidewire is a stainless-steel guidewire with a fixed distal-spring coil and is designed to deliver exceptional support for therapeutic devices.



Features

Deliverability

Fixed distal silicone-coated spring coil eases the crossing of difficult lesions.

Control

Core-to-tip design and medium-stiff coil facilitates superior torque transmission and control.

Support

Stainless-steel mandrel (0.014 inch) is designed to provide excellent support.


Cardio Flow

800.294.5517 ▪ cardioflow.net ▪ customerservice@cardioflow.net

Mailing address: 525 Main Street, Box 120018 ▪ St. Paul, MN 55112 USA

©2022 Cardio Flow, Inc. All rights reserved. 8032 Rev B

FreedomFlow® Peripheral Guidewire



Cardio Flow FreedomFlow Peripheral Guidewire GW1001 Specifications

Model	Nominal Core Outer Diameter	Maximum Core Outer Diameter	Maximum Spring Tip Outer Diameter	Recommended Introducer Sheath	Length	Tip Load Force
GW1001	0.013 inch (0.33 mm)	0.014 inch (0.36 mm)	0.0144 inch (0.366 mm)	0.014 inch minimum (0.36 mm) lumen compatibility	325 cm	13 gf (Reference)



FreedomFlow®
by Cardio Flow
PERIPHERAL GUIDEWIRE

To order

Please contact Cardio Flow at 800-294-5517 to order the FreedomFlow Peripheral Guidewire.

Important information

CAUTION: Federal (USA) law restricts this device to sale by or on the order of a physician.

INDICATIONS FOR USE: The Cardio Flow FreedomFlow Peripheral Guidewire is intended for temporary placement in peripheral vasculature to facilitate the placement and exchange of diagnostic and therapeutic devices during percutaneous intravascular procedures. This guidewire device is intended for peripheral vascular use only.

Refer to the Instructions for Use for a complete listing of the indications, contraindications, warnings, and precautions.

Cardio Flow

800.294.5517 ▪ cardioflow.net ▪ customerservice@cardioflow.net

Mailing address: 525 Main Street, Box 120018 ▪ St. Paul, MN 55112 USA

©2022 Cardio Flow, Inc. All rights reserved. 8032 Rev B